

BMP-3 Polyclonal Antibody

Catalog No: #40649

Package Size: #40649-1 50ul #40649-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	BMP-3 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC
Species Reactivity	Hu Ms Rt
Specificity	BMP-3 Polyclonal Antibody detects endogenous levels of BMP-3 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human BMP-3.
Target Name	BMP-3
Other Names	BMP3; BMP3A; Bone morphogenetic protein 3; BMP-3; Bone morphogenetic protein 3A; BMP-3A; Osteogenin
Accession No.	Swiss-Prot: P12645NCBI Gene ID: 651
Uniprot	P12645
GeneID	651;
SDS-PAGE MW	53kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

Images



Western Blot analysis of Jurkat cells using BMP-3 Polyclonal Antibody

Note: This product is for in vitro research use only