

Transthyretin

Catalog No: #AP78636

Package Size: #AP78636-1 50ug #AP78636-2 100ug #AP78636-3 1mg


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Transthyretin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	21-147AA
Other Names	PALB,ATTR,Prealbumin,TBPA
Accession No.	P02766Gene name:TTR
Uniprot	P02766
GeneID	7276;
Calculated MW	13.97
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Images



Background

Thyroid hormone-binding protein. Probably transports thyroxine from the bloodstream to the brain.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only