

BMP-5 Polyclonal Antibody

Catalog No: #40650

Package Size: #40650-1 50ul #40650-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	BMP-5 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms
Specificity	BMP-5 Polyclonal Antibody detects endogenous levels of BMP-5 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human BMP-5.
Target Name	BMP-5
Other Names	BMP5; Bone morphogenetic protein 5; BMP-5
Accession No.	Swiss-Prot: P22003NCBI Gene ID: 653
Uniprot	P22003
GeneID	653;
SDS-PAGE MW	52kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.


ELISA: 1/40000.

Not yet tested in other applications.

Images



Western Blot analysis of HeLa cells using BMP-5 Polyclonal Antibody



Western Blot analysis of VEC 293 A549 cells using BMP-5 Polyclonal Antibody

Note: This product is for in vitro research use only