

Complement receptor type 2

Catalog No: #AP78644



Package Size: #AP78644-1 50ug #AP78644-2 100ug #AP78644-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Complement receptor type 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Mouse
Immunogen Description	517-587AA
Other Names	Complement C3d receptor,CD21
Accession No.	P19070Gene name:Cr2
Uniprot	P19070
Calculated MW	7.81
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Receptor for complement C3d. Participates in B lymphocytes activation.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only