

Hepatitis A virus cellular receptor 1 homolog

Catalog No: #AP78659



Package Size: #AP78659-1 50ug #AP78659-2 100ug #AP78659-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Hepatitis A virus cellular receptor 1 homolog
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rat
Immunogen Description	22-235AA
Other Names	Kim1,Kidney injury molecule 1,T cell immunoglobulin and mucin domain-containing protein 1,KIM-1,TIMD-1
Accession No.	O54947Gene name:Havcr1
Uniprot	O54947
GeneID	286934;
Calculated MW	23.54
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

May play a role in T-helper cell development and the regulation of asthma and allergic diseases. Receptor for TIMD4. May play a role in kidney injury and repair.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only