

BRCA1 Polyclonal Antibody

Catalog No: #40655

Package Size: #40655-1 50ul #40655-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	BRCA1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF ELISA
Species Reactivity	Hu Rt
Specificity	BRCA1 Polyclonal Antibody detects endogenous levels of BRCA1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human BRCA1 around the non-phosphorylation site of S1423.
Target Name	BRCA1
Other Names	BRCA1; RNF53; Breast cancer type 1 susceptibility protein; RING finger protein 53
Accession No.	Swiss-Prot: P38398NCBI Gene ID: 672
Uniprot	P38398
GeneID	672;
SDS-PAGE MW	207kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000.


ELISA: 1/5000.

Not yet tested in other applications.

Images



Western Blot analysis of HT29 cells using BRCA1 Polyclonal Antibody



Western Blot analysis of 293T-UV 22RV1 cells using BRCA1 Polyclonal Antibody

Note: This product is for in vitro research use only