

# Prolactin

Catalog No: #AP78693



Package Size: #AP78693-1 50ug #AP78693-2 100ug #AP78693-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Prolactin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Immunogen Description	31-229AA
Other Names	null
Accession No.	P01239Gene name:PRL
Uniprot	P01239
GeneID	280901;
Calculated MW	21.89
Tag Info	His
Formulation	50mM NaH <sub>2</sub> PO <sub>4</sub> , 500mM NaCl Buffer with 500mM Imidazole, 10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

## Background

Prolactin acts primarily on the mammary gland by promoting lactation.

## References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only