

# Resistin

Catalog No: #AP78694

Package Size: #AP78694-1 50ug #AP78694-2 100ug #AP78694-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Resistin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	Recombinant protein
Other Names	UNQ407/PRO1199,FIZZ3,HXCP1,RSTN,Adipose tissue-specific secretory factor,C/EBP-epsilon-regulated myeloid-specific secreted cysteine-rich protein,Cysteine-rich secreted protein A12-alpha-like 2,Cysteine-rich secreted protein FIZZ3,ADSF
Accession No.	Q9HD89Gene name:RETN
Uniprot	Q9HD89
GeneID	56729;
Tag Info	His
Formulation	50mM NaH <sub>2</sub> PO <sub>4</sub> , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

## Background

Hormone that seems to suppress insulin ability to stimulate glucose uptake into adipose cells. Potentially links obesity to diabetes.

## References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only