

## Resistin

Catalog No: #AP78695



Package Size: #AP78695-1 50ug #AP78695-2 100ug #AP78695-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Resistin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Mouse
Immunogen Description	21-114aa
Other Names	Fizz3, Adipose tissue-specific secretory factor, Adipose-specific cysteine-rich secreted protein A12-alpha, Cysteine-rich secreted protein FIZZ3, ADSF
Accession No.	Q99P87 Gene name: Retn
Uniprot	Q99P87
GeneID	57264;
Calculated MW	10.34
Tag Info	His
Formulation	50mM NaH <sub>2</sub> PO <sub>4</sub> , 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.) Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

## Background

Hormone that seems to suppress insulin ability to stimulate glucose uptake into adipose cells. Potentially links obesity to diabetes.

## References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only