

Coagulation factor XI

Catalog No: #AP78696



Package Size: #AP78696-1 50ug #AP78696-2 100ug #AP78696-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Coagulation factor XI
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	388-625AA
Other Names	Plasma thromboplastin antecedent,PTA
Accession No.	P03951Gene name:F11
Uniprot	P03951
GeneID	2160;
Calculated MW	26.18
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Factor XI triggers the middle phase of the intrinsic pathway of blood coagulation by activating factor IX.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only