

Uteroglobin

Catalog No: #AP78698



Package Size: #AP78698-1 50ug #AP78698-2 100ug #AP78698-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Uteroglobin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	CEICPSFQRVIETLLMDTPSSYEAA
Other Names	CC10,CCSP,UGB,Clara cell phospholipid-binding protein,Clara cells 10 kDa secretory protein,Secretoglobin family 1A member 1,Urinary protein 1,CCPBP,CC10,UP-1
Accession No.	P11684Gene name:SCGB1A1
Uniprot	P11684
GeneID	7356;
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Binds phosphatidylcholine, phosphatidylinositol, polychlorinated biphenyls (PCB) and weakly progesterone, potent inhibitor of phospholipase A2.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only