

Pro-neuropeptide Y

Catalog No: #AP78709



Package Size: #AP78709-1 50ug #AP78709-2 100ug #AP78709-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Pro-neuropeptide Y
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	29-64AA
Other Names	null
Accession No.	P01303Gene name:NPY
Uniprot	P01303
GeneID	4852;
Calculated MW	3.96
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole, 10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

NPY is implicated in the control of feeding and in secretion of gonadotrophin-release hormone.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only