

Gastric inhibitory polypeptide

Catalog No: #AP78710



Package Size: #AP78710-1 50ug #AP78710-2 100ug #AP78710-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------------|--|
| Product Name | Gastric inhibitory polypeptide |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Human |
| Immunogen Description | CAQKGKKNDWKHNITQ |
| Other Names | Glucose-dependent insulinotropic polypeptide, Incretin hormone |
| Accession No. | P09681 Gene name: GIP |
| Uniprot | P09681 |
| GeneID | 2695; |
| Tag Info | His |
| Formulation | 50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) |
| Storage | Store at -20C. (Avoid repeated freezing and thawing.) Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |

Background

Potent stimulator of insulin secretion and relatively poor inhibitor of gastric acid secretion.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only