

Serum amyloid A protein

Catalog No: #AP78712



Package Size: #AP78712-1 50ug #AP78712-2 100ug #AP78712-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Serum amyloid A protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Immunogen Description	19-130AA
Other Names	SAA
Accession No.	P35541Gene name:SAA1
Uniprot	P35541
GeneID	506412;
Calculated MW	12.32
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole, 10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Major acute phase reactant. Apolipoprotein of the HDL complex.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only