

## Serum amyloid A-1 protein

Catalog No: #AP78713


Package Size: #AP78713-1 50ug #AP78713-2 100ug #AP78713-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Serum amyloid A-1 protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	19-122aa
Other Names	null
Accession No.	P0DJ18Gene name:SAA1
Uniprot	P0DJ18
GeneID	6288;
Calculated MW	11.44
Tag Info	His
Formulation	50mM NaH <sub>2</sub> PO <sub>4</sub> , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

## Images



## Background

Major acute phase protein.

## References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

---

Note: This product is for in vitro research use only