

Amyloid beta A4 protein

Catalog No: #AP78722



Package Size: #AP78722-1 50ug #AP78722-2 100ug #AP78722-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Amyloid beta A4 protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	CNAQFRHNSGYQVHHQKLVFFAQNVGSNKGAIIGLMVGGVV
Other Names	A4,AD1,ABPP,APPI,Alzheimer disease amyloid protein,Amyloid precursor protein,Beta-amyloid precursor protein,Cerebral vascular amyloid peptide,PreA4,Protease nexin-II,APP,CVAP,PN-II
Accession No.	P05067Gene name:APP
Uniprot	P05067
GeneID	351;
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

N-APP binds TNFRSF21 triggering caspase activation and degeneration of both neuronal cell bodies (via caspase-3) and axons (via caspase-6).

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only