

C/EBP β Polyclonal Antibody

Catalog No: #40658

Package Size: #40658-1 50ul #40658-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	C/EBP β Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms
Specificity	C/EBP β Polyclonal Antibody detects endogenous levels of C/EBP β protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human C/EBP β .
Target Name	C/EBP β
Other Names	CEBPB; LAP; TCF5; PP9092; CCAAT/enhancer-binding protein beta; C/EBP beta; Liver activator protein; Nuclear factor NF-IL6; Transcription factor 5; TCF-5
Accession No.	Swiss-Prot: P17676NCBI Gene ID: 1051
Uniprot	P17676
GeneID	1051;
SDS-PAGE MW	49kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/5000.

Not yet tested in other applications.

Images

Western Blot analysis of 293T cells using C/EBP β Polyclonal Antibody



Note: This product is for in vitro research use only