

## Glypican-3

Catalog No: #AP78726

Package Size: #AP78726-1 50ug #AP78726-2 100ug #AP78726-3 1mg


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Glypican-3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	25-206AA
Other Names	OCI5,GTR2-2,Intestinal protein OCI-5,MXR7
Accession No.	P51654Gene name:GPC3
Uniprot	P51654
GeneID	2719;
Calculated MW	20.02
Tag Info	His
Formulation	50mM NaH <sub>2</sub> PO <sub>4</sub> , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

## Images



## Background

Cell surface proteoglycan that bears heparan sulfate. Inhibits the dipeptidyl peptidase activity of DPP4. May be involved in the suppression/modulation of growth in the predominantly mesodermal tissues and organs. May play a role in the modulation of IGF2 interactions with its receptor and thereby modulate its function. May regulate growth and tumor predisposition.

## References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

---

Note: This product is for in vitro research use only