

T-cell surface glycoprotein CD3 gamma chain

Catalog No: #AP78732

Package Size: #AP78732-1 50ug #AP78732-2 100ug #AP78732-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	T-cell surface glycoprotein CD3 gamma chain
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	TEDKKKWNLG SNAKDPRGMYQC
Other Names	T3G, T-cell receptor T3 gamma chain, CD3g
Accession No.	P09693 Gene name: CD3G
Uniprot	P09693
GeneID	917;
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.) Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

The CD3 complex mediates signal transduction.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only