

# Leukotriene C4 synthase

Catalog No: #AP78737



Package Size: #AP78737-1 50ug #AP78737-2 100ug #AP78737-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Leukotriene C4 synthase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	1-148AA
Other Names	Leukotriene-C(4) synthase
Accession No.	Q16873Gene name:LTC4S
Uniprot	Q16873
GeneID	4056;
Calculated MW	16.28
Tag Info	His
Formulation	50mM NaH <sub>2</sub> PO <sub>4</sub> , 500mM NaCl Buffer with 500mM Imidazole, 10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

## Background

Catalyzes the conjugation of leukotriene A4 with reduced glutathione to form leukotriene C4.

## References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only