

C1q-B Polyclonal Antibody

Catalog No: #40660

Package Size: #40660-1 50ul #40660-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	C1q-B Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF ELISA
Species Reactivity	Hu
Specificity	C1q-B Polyclonal Antibody detects endogenous levels of C1q-B protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human C1q-B.
Target Name	C1q-B
Other Names	C1QB; Complement C1q subcomponent subunit B
Accession No.	Swiss-Prot: P02746NCBI Gene ID: 713
Uniprot	P02746
GeneID	713;
SDS-PAGE MW	28kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.


Immunofluorescence: 1/200 - 1/1000.

ELISA: 1/40000.

Not yet tested in other applications.

Images

Western Blot analysis of Jurkat cells using C1q-B Polyclonal Antibody



Note: This product is for in vitro research use only