

6-phosphofructokinase, liver type

Catalog No: #AP78746



Package Size: #AP78746-1 50ug #AP78746-2 100ug #AP78746-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	6-phosphofructokinase, liver type
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rat
Immunogen Description	401-780AA
Other Names	Pfk-I, Phosphofructo-1-kinase isozyme B, Phosphofructokinase 1, Phosphohexokinase, PFK-B
Accession No.	P30835 Gene name: PfkI
Uniprot	P30835
GeneID	25741;
Calculated MW	41.8
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.) Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

null

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only