

Beta-lactoglobulin

Catalog No: #AP78749



Package Size: #AP78749-1 50ug #AP78749-2 100ug #AP78749-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Beta-lactoglobulin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Immunogen Description	Recombinant protein
Other Names	null
Accession No.	P02754Gene name:LGB
Uniprot	P02754
GeneID	280838;
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Primary component of whey, it binds retinol and is probably involved in the transport of that molecule.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only