

## Beclin-1

Catalog No: #AP78751



Package Size: #AP78751-1 50ug #AP78751-2 100ug #AP78751-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                       |   |
|-----------------------|---|
| Product Name          | Beclin-1  |
| Brief Description     | Recombinant Protein   |
| Host Species          | E.coli  |
| Purification          | Greater than 90% by SDS-PAGE  |
| Species Reactivity    | Human   |
| Immunogen Description | 150-450AA   |
| Other Names           | GT197,Coiled-coil myosin-like BCL2-interacting protein,Protein GT197  |
| Accession No.         | Q14457Gene name:BECN1   |
| Uniprot               | Q14457  |
| GeneID                | 8678;   |
| Calculated MW         | 33.11   |
| Tag Info              | His   |
| Formulation           | 50mM NaH <sub>2</sub> PO <sub>4</sub> , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)   |
| Storage               | Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended.<br>Store working aliquots at 4°C for up to one week. |

## Background

Beclin-1-C 35 kDa localized to mitochondria can promote apoptosis; it induces the mitochondrial translocation of BAX and the release of proapoptotic factors.

## References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only