

Growth hormone receptor

Catalog No: #AP78763



Package Size: #AP78763-1 50ug #AP78763-2 100ug #AP78763-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Growth hormone receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	CGYVSTDQLNKIMP
Other Names	Somatotropin receptor
Accession No.	P10912Gene name:GHR
Uniprot	P10912
GeneID	2690;
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Isoform 2 up-regulates the production of GHBP and acts as a negative inhibitor of GH signaling.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only