

## Vimentin

Catalog No: #AP78769



Package Size: #AP78769-1 50ug #AP78769-2 100ug #AP78769-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Vimentin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	269-407AA
Other Names	null
Accession No.	P08670Gene name:VIM
Uniprot	P08670
GeneID	7431;
Calculated MW	15.29
Tag Info	His
Formulation	50mM NaH <sub>2</sub> PO <sub>4</sub> , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

## Background

Involved with LARP6 in the stabilization of type I collagen mRNAs for CO1A1 and CO1A2.

## References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only