

# Bile salt export pump

Catalog No: #AP78770

Package Size: #AP78770-1 50ug #AP78770-2 100ug #AP78770-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Bile salt export pump
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Mouse
Immunogen Description	1033-1321AA
Other Names	Bsep,Spgp,ATP-binding cassette sub-family B member 11,Sister of P-glycoprotein
Accession No.	Q9QY30Gene name:Abcb11
Uniprot	Q9QY30
GeneID	27413;
Calculated MW	31.79
Tag Info	His
Formulation	50mM NaH <sub>2</sub> PO <sub>4</sub> , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

## Background

Involved in the ATP-dependent secretion of bile salts into the canalculus of hepatocytes.

## References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only