

## Abl1 Polyclonal Antibody

Catalog No: #40663

Package Size: #40663-1 50ul #40663-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Abl1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt Mk
Specificity	Abl1 Polyclonal Antibody detects endogenous levels of Abl1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human Abl1 around the non-phosphorylation site of T735.
Target Name	Abl1
Other Names	ABL1; ABL; JTK7; Tyrosine-protein kinase ABL1; Abelson murine leukemia viral oncogene homolog 1; Abelson tyrosine-protein kinase 1; Proto-oncogene c-Abl; p150
Accession No.	Swiss-Prot: P00519NCBI Gene ID: 25
Uniprot	P00519
GeneID	25;
SDS-PAGE MW	140kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


## Application Details

Western Blot: 1/500 - 1/2000.


ELISA: 1/20000.

Not yet tested in other applications.

## Images



Western Blot analysis of COS cells using Abl1 Polyclonal Antibody



Western Blot analysis of cos-7 cells using AbI1 Polyclonal Antibody

Note: This product is for in vitro research use only