

Interferon-induced 35 kDa protein homolog

Catalog No: #AP78775

Package Size: #AP78775-1 50ug #AP78775-2 100ug #AP78775-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Interferon-induced 35 kDa protein homolog
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Mouse
Immunogen Description	2-286AA
Other Names	null
Accession No.	Q9D8C4Gene name:lfi35
Uniprot	Q9D8C4
GeneID	70110;
Calculated MW	31.35
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

null

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only