

Cathelicidin antimicrobial peptide

Catalog No: #AP78776



Package Size: #AP78776-1 50ug #AP78776-2 100ug #AP78776-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cathelicidin antimicrobial peptide
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	31-131AA
Other Names	HSD26,CAP18,FALL39,18 kDa cationic antimicrobial protein,CAP-18
Accession No.	P49913Gene name:CAMP
Uniprot	P49913
GeneID	820;
Calculated MW	11.11
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Binds to bacterial lipopolysaccharides (LPS), has antibacterial activity.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only