

Serine protease inhibitor A3N

Catalog No: #AP78777



Package Size: #AP78777-1 50ug #AP78777-2 100ug #AP78777-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Serine protease inhibitor A3N
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Mouse
Immunogen Description	21-418AA
Other Names	Spi2
Accession No.	Q91WP6Gene name:Serpina3n
Uniprot	Q91WP6
Calculated MW	43.78
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

null

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only