

Tryptase gamma

Catalog No: #AP78778



Package Size: #AP78778-1 50ug #AP78778-2 100ug #AP78778-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tryptase gamma
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	38-321AA
Other Names	PRSS31, TMT, Serine protease 31, Transmembrane trypsin
Accession No.	Q9NRR2 Gene name:TPSG1
Uniprot	Q9NRR2
GeneID	25823;
Calculated MW	31.24
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.) Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

null

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only