

## Biglycan

Catalog No: #AP78789



Package Size: #AP78789-1 50ug #AP78789-2 100ug #AP78789-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Biglycan
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Immunogen Description	38-369AA
Other Names	Bone/cartilage proteoglycan I,PG-S1
Accession No.	O02678Gene name:BGN
Uniprot	O02678
GeneID	403905;
Calculated MW	36.52
Tag Info	His
Formulation	50mM NaH <sub>2</sub> PO <sub>4</sub> , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

## Background

May be involved in collagen fiber assembly.

## References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only