

Cadherin-19 Polyclonal Antibody

Catalog No: #40666

Package Size: #40666-1 50ul #40666-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cadherin-19 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF
Species Reactivity	Hu
Specificity	Cadherin-19 Polyclonal Antibody detects endogenous levels of Cadherin-19 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human Cadherin-19.
Target Name	Cadherin-19
Other Names	CDH19; CDH7L2; Cadherin-19
Accession No.	Swiss-Prot: Q9H159NCBI Gene ID: 28513
Uniprot	Q9H159
GeneID	28513;
SDS-PAGE MW	87kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000.

Images



Western Blot analysis of COLO205 cells using Cadherin-19 Polyclonal Antibody

Note: This product is for in vitro research use only