

Interferon alpha-inducible protein 6

Catalog No: #AP78814



Package Size: #AP78814-1 50ug #AP78814-2 100ug #AP78814-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------------|--|
| Product Name | Interferon alpha-inducible protein 6 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Human |
| Immunogen Description | Recombinant protein |
| Other Names | G1P3, Interferon-induced protein 6-16, Irfi-6-16 |
| Accession No. | P09912 Gene name: IFI6 |
| Uniprot | P09912 |
| GeneID | 2537; |
| Tag Info | His |
| Formulation | 50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (PH8.0) |
| Storage | Store at -20C. (Avoid repeated freezing and thawing.) Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |

Background

null

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only