

## Arginase-1

Catalog No: #AP78819


Package Size: #AP78819-1 50ug #AP78819-2 100ug #AP78819-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Arginase-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Mouse
Immunogen Description	1-323aa
Other Names	Liver-type arginase, Type I arginase
Accession No.	Q61176 Gene name: Arg1
Uniprot	Q61176
GeneID	11846;
Calculated MW	35.53
Tag Info	His
Formulation	50mM NaH <sub>2</sub> PO <sub>4</sub> , 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.) Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

## Images



## Background

null

## References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

---

Note: This product is for in vitro research use only