

Atrial natriuretic peptide receptor 2

Catalog No: #AP78827



Package Size: #AP78827-1 50ug #AP78827-2 100ug #AP78827-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Atrial natriuretic peptide receptor 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	SLRRSSCFGGRMDRIGASGLGCNSFRY
Other Names	ANPRB,Atrial natriuretic peptide receptor type B,Guanylate cyclase B,ANP-B,GC-B
Accession No.	P20594Gene name:NPR2
Uniprot	P20594
GeneID	4882;
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Receptor for the C-type natriuretic peptide NPPC/CNP hormone. Has guanylate cyclase activity upon binding of its ligand. May play a role in the regulation of skeletal growth.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only