

Bone morphogenetic protein 3B

Catalog No: #AP78835



Package Size: #AP78835-1 50ug #AP78835-2 100ug #AP78835-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Bone morphogenetic protein 3B
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	369-478aa
Other Names	BMP3B, Bone-inducing protein, Growth/differentiation factor 10, BIP, GDF-10
Accession No.	P55107 Gene name: GDF10
Uniprot	P55107
GeneID	2662;
Calculated MW	12.1
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0)
Storage	Store at -20°C. (Avoid repeated freezing and thawing.) Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

null

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only