

Huntingtin

Catalog No: #AP78844



Package Size: #AP78844-1 50ug #AP78844-2 100ug #AP78844-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Huntingtin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	1-115AA
Other Names	HD,IT15,Huntington disease protein,HD protein
Accession No.	P42858Gene name:HTT
Uniprot	P42858
GeneID	3064;
Calculated MW	12.65
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

May play a role in microtubule-mediated transport or vesicle function.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only