

Lutropin subunit beta

Catalog No: #AP78848



Package Size: #AP78848-1 50ug #AP78848-2 100ug #AP78848-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Lutropin subunit beta
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Immunogen Description	21-141AA
Other Names	Luteinizing hormone subunit beta,LH-B
Accession No.	P01232Gene name:LHB
Uniprot	P01232
GeneID	397153;
Calculated MW	13.31
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Promotes spermatogenesis and ovulation by stimulating the testes and ovaries to synthesize steroids.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only