

Gamma-crystallin D

Catalog No: #AP78849



Package Size: #AP78849-1 50ug #AP78849-2 100ug #AP78849-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Gamma-crystallin D
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	Recombinant protein
Other Names	CRYG4, Gamma-D-crystallin, Gamma-crystallin 4
Accession No.	P07320 Gene name: CRYGD
Uniprot	P07320
GeneID	1421;
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.) Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Crystallins are the dominant structural components of the vertebrate eye lens.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only