

Beta-Ala-His dipeptidase

Catalog No: #AP78854



Package Size: #AP78854-1 50ug #AP78854-2 100ug #AP78854-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Beta-Ala-His dipeptidase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	Recombinant Protein
Other Names	UNQ1915/PRO4380,CN1,CPGL2,CNDP dipeptidase 1,Carnosine dipeptidase 1,Glutamate carboxypeptidase-like protein 2,Serum carnosinase
Accession No.	Q96KN2Gene name:CNDP1
Uniprot	Q96KN2
GeneID	84735;
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

null

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only