

Metal regulatory transcription factor 1

Catalog No: #AP78862



Package Size: #AP78862-1 50ug #AP78862-2 100ug #AP78862-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Metal regulatory transcription factor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	588-707aa
Other Names	MRE-binding transcription factor, Transcription factor MTF-1
Accession No.	Q14872 Gene name: MTF1
Uniprot	Q14872
GeneID	4520;
Calculated MW	13.2
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.) Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Activates the metallothionein I promoter. Binds to the metal responsive element (MRE).

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only