

## Glycylpeptide N-tetradecanoyltransferase 1

Catalog No: #AP78865



Package Size: #AP78865-1 50ug #AP78865-2 100ug #AP78865-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Glycylpeptide N-tetradecanoyltransferase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	1-80AA
Other Names	NMT, Myristoyl-CoA:protein N-myristoyltransferase 1, Peptide N-myristoyltransferase 1, NMT 1
Accession No.	P30419 Gene name: NMT1
Uniprot	P30419
GeneID	4836;
Calculated MW	8.8
Tag Info	His
Formulation	50mM NaH <sub>2</sub> PO <sub>4</sub> , 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.) Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

## Background

Adds a myristoyl group to the N-terminal glycine residue of certain cellular and viral proteins.

## References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only