

Transforming growth factor beta-3

Catalog No: #AP78874



Package Size: #AP78874-1 50ug #AP78874-2 100ug #AP78874-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------------|---|
| Product Name | Transforming growth factor beta-3 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Human |
| Immunogen Description | Recombinant protein |
| Other Names | null |
| Accession No. | P10600Gene name:TGFB3 |
| Uniprot | P10600 |
| GeneID | 7043; |
| Tag Info | His |
| Formulation | 50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0) |
| Storage | Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |

Background

Involved in embryogenesis and cell differentiation.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only