

Hyaluronidase-1

Catalog No: #AP78876

Package Size: #AP78876-1 50ug #AP78876-2 100ug #AP78876-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Hyaluronidase-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	183-435aa
Other Names	LUCA1,Hyaluronoglucosaminidase-1,Lung carcinoma protein 1,LuCa-1
Accession No.	Q12794Gene name:HYAL1
Uniprot	Q12794
GeneID	3373;
Calculated MW	27.83
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

May have a role in promoting tumor progression. May block the TGFB1-enhanced cell growth.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only