

## Stanniocalcin-2

Catalog No: #AP78878

Package Size: #AP78878-1 50ug #AP78878-2 100ug #AP78878-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Stanniocalcin-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	25-302AA
Other Names	Stanniocalcin-related protein,STC-related protein
Accession No.	O76061Gene name:STC2
Uniprot	O76061
GeneID	8614;
Calculated MW	30.58
Tag Info	His
Formulation	50mM NaH <sub>2</sub> PO <sub>4</sub> , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

## Background

Has an anti-hypocalcemic action on calcium and phosphate homeostasis.

## References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only