

Neural cell adhesion molecule 1

Catalog No: #AP78882



Package Size: #AP78882-1 50ug #AP78882-2 100ug #AP78882-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Neural cell adhesion molecule 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	CTKDGEQIEQEEDDEKYIFSDDSSQL
Other Names	NCAM,CD56
Accession No.	P13591Gene name:NCAM1
Uniprot	P13591
GeneID	4684;
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

(Microbial infection) Acts as a receptor for rabies virus.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only