

Glutamate decarboxylase 1

Catalog No: #AP78908



Package Size: #AP78908-1 50ug #AP78908-2 100ug #AP78908-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Glutamate decarboxylase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	Recombinant Protein
Other Names	GAD,GAD67,67 kDa glutamic acid decarboxylase,Glutamate decarboxylase 67 kDa isoform,GAD-67
Accession No.	Q99259Gene name:GAD1
Uniprot	Q99259
GeneID	2571;
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Catalyzes the production of GABA.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only