

## Humanin-like 6

Catalog No: #AP78909

Package Size: #AP78909-1 50ug #AP78909-2 100ug #AP78909-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Humanin-like 6
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	1-24AA
Other Names	MT-RNR2-like protein 6
Accession No.	P0CJ73Gene name:MTRNR2L6
Uniprot	P0CJ73
GeneID	100463482;
Calculated MW	2.64
Tag Info	His
Formulation	50mM NaH <sub>2</sub> PO <sub>4</sub> , 500mM NaCl Buffer with 500mM Imidazole, 10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

## Background

Plays a role as a neuroprotective and antiapoptotic factor.

## References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only